

# NYLON TAPE

## - STERILE -

389907R01



01/2014  
LAB100144760V1



0086

ar	NYLON TAPE شريط نايلون	et	NAILONLING	ko	나일론 테이프	ru	НЕЙЛОНОВАЯ ЛЕНТА
bg	НАЙЛОНОВА ЛЕНТА	fi	NAILONNAUHA	lt	NAILONO JUOSTELĖ	sk	NYLONOVÁ PÁSKA
cs	NYLONOVÁ PÁSKA	fr	BANDE NYLON	lv	NEILONA LENTE	sl	NAJLONSKI TRAK
da	NYLONTAPE	hr	NAJLONSKA VRPČA	nl	NYLON BAND	sv	NYLONBAND
de	NYLONBAND	hu	NYLON SZALAG	no	NYLONTAPE	tr	NAYLON BANT
el	ΝΑΙΛΟΝ ΤΑΙΝΙΑ	id	PLESTER NILON	pl	TAŚMA NYLONOWA	zh-cn	尼龙带
en	NYLON TAPE	it	NASTRO DI NYLON	pt	FITA NYLON	zh-tw	尼龍膠帶
es	CINTA DE NYLON	kk	НЕЙЛОН БАУЫ	ro	BANDĂ DE NYLON		

### Instructions for use

## NYLON TAPE

### STERILE

#### DESCRIPTION

Nylon Tape is a white (undyed) flat woven or braided Nylon-(Polyamide 6,6) tape. Nylon Tape is available in various widths and lengths, non-needed only.

#### PERFORMANCE

Because of its flat woven structure, Nylon Tape is indicated for temporary ligation or retraction of tissues, organs or other anatomical structures during surgical procedures without injuring the tissue involved.

#### INDICATIONS

Nylon Tape is indicated for temporary ligation or retraction of tissues, organs or other anatomical structures during surgical procedures. The device is not intended to be implanted in the body.

#### APPLICATION

Nylon Tape is used for temporary retraction and/or fixing of organs or parts of organs (e.g. intestine) during an operation. After completion of the operation the tape should be removed and disposed of.

#### CONTRAINDICATIONS

There are no known contraindications. However, Nylon Tape is not indicated for implantation.

#### WARNINGS / PRECAUTIONS

Care should be taken during tissue retraction to avoid injury. Nylon Tape is provided as a sterile product. Do not resterilize. Do not use if packaging is opened or damaged. Discard opened unused products.

Reuse of this device (or portions of this device) may create a risk of product degradation and cross-contamination, which may lead to failure of treatment, infection or transmission of bloodborne pathogens to patients and users.

Safety and effectiveness of Nylon Tape have not been established for use as an implanted device. The risks associated with product implantation have not been established.

#### ADVERSE REACTIONS

None known.

#### STERILITY

Nylon Tape is sterilized by irradiation. Do not resterilize!  
Do not use if package is opened or damaged! Discard open, unused product.

#### STORAGE

No special storage conditions required. Do not use after expiry date!

#### SYMBOLS USED ON LABELING

	Sterilized using irradiation		Do not reuse
	Do not use if package is damaged		Do not resterilize
	Batch code		Use by - year and month
	CE-mark and Identification Number of Notified Body. The product meets the essential requirements of Medical Device Directive 93/42/EEC		Manufacturer
	Authorized representative in the European Community		Caution: See instructions for use
			Number of Sachets in Box
			Catalogue Number



Ethicon, LLC  
475 C Street  
Los Frailes Industrial Park  
Suite 401  
Guaynabo, PR 00969 USA  
+1-513-337-6928



Johnson & Johnson Medical Ltd  
Simpson Parkway  
Kirkton Campus  
Livingston  
EH54 7AT  
United Kingdom

© Ethicon, Inc. 2014