Rocket

Verres Insufflation Needles

INSTRUCTIONS FOR USE

Scope:

These instructions cover all Rocket Verres pneumoperitoneal insufflation needles. R57920, R57921 and derivatives

Figure 1: Rocket Verres Needle

Indications: For trans-umbilical insufflation of CO₂ or similar gas media into the abdominal cavity prior to endoscopy. This device should only be used by or under the supervision of trained personnel and in conjunction with current clinical practice guidelines such as those published by the Royal College of Surgeons, the Royal College of Obstetricians & Gynaecologists and NICE.

Contraindications: Not intended for use where endoscopic techniques are contraindicated such as those where abdominal architecture is uncertain or cannot be established by ultrasound or similar imaging techniques.

Preparation:

1. Follow normal local hospital policy for preparation for laparoscopy; such as the routine cleansing of the skin and perineal area with a suitable antiseptic solution.
2. The patient's bladder should be drained prior to insertion of the needle.
3. The patient may be placed in a modified lithotomy position with the legs flexed to 45 degrees and with 15 degrees of Trendelenberg (head down) tilt.

Pre-operative:

4. Prior to insertion, attach a suitable tubing set with filter such as the Rocket Filter set to the insufflation source.
5. Attach the tubing set to the luer lock of the needle stopcock.

NOTE: The orientation of the stopcock - the tap is *open* (on) when the arms are in line with the needle.

6. Pass 1-2 litres of gas through the needle, there should be no appreciable obstruction to flow.


Figure 2: Rocket Filter Set

Recommended Operative Technique:


1. Make a 1-2cm incision is just below the umbilical margin, insert the needle through the incision whilst using the other hand to lift the abdominal wall away from the underlying viscera and blood vessels (). Insertion is ideally made through the deepest part of the umbilicus as the thickness of the abdominal wall is minimal and the peritoneum is closely adherent.

2. During insertion through the abdominal wall, the hub window will show RED to indicate the needle tip is exposed.

3. As the needle passes through the peritoneum a distinct 'click' can be heard as the spring loaded obturator snaps forward to protect the internal organs from the needle point.

4. Check the hub window is showing GREEN to indicate the obturator is fully forward and protecting the needle tip.

5. A RED indicator showing or partially showing, indicates that the obturator is retracted and the needle tip may be exposed. In this condition, remove the needle, ensure that there is no tissue obstructing the free movement of the obturator and repeat the insertion procedure.

Aspiration Test:

1. Inject 10ml Normal saline through the insufflation needle into the cavity. Attempt to withdraw the fluid. If the needle lies correctly in the peritoneal cavity, no fluid should be aspirated.

2. If clear fluid is freely aspirated, then the needle may still lie in the abdominal wall and has not punctured the peritoneum. The presence of blood or bowel contents in the aspirated fluid indicates that the needle has punctured a blood vessel or loop of bowel and immediate exploratory and corrective action should be taken.

3. Inflation of the peritoneal cavity with between 1000 - 4000ml of CO2, nitrous oxide or occasionally nitrogen at 15-20mmHg from the Gas Insufflator may follow once the correct position is confirmed.

CAUTION: The Rocket Verres Needle should only be used for primary insufflation

Disposal

This device should be handled and disposed of in accordance with local hospital policy and with regard to all applicable regulations, including but without limitation to, 'those pertaining to human health & safety and care of the environment.

This device is not manufactured with natural rubber latex

ROCKET MEDICAL PLC Sedling Road, Washington, England, NE38 9BZ www.rocketmedical.com

Rocket Medical GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

DO NOT RESTERILISE

Unless opened or damaged, contents of package are sterile. Sterile unless package has been opened or damaged

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Copyright© 1998-2020