



Push Button

Blood Collection Set




- REF 367335
- REF 367336
- REF 367338
- REF 367341
- REF 367342
- REF 367344

- REF 367323
- REF 367324
- REF 367326

For Blood Collection

Note: Should be performed wearing gloves.


1. Thread the Luer Adapter into Holder.
2. Remove needle sheath.
3. (a) Insert needle in vein by holding wings. (b) Observe for the presence of blood in the chamber. Collect the blood specimen according to your facility's procedure.
4. Retract the Needle: Depress the button. The needle will slide out of the venipuncture site and lock into place. Do not impede device retraction. Cover the puncture site with a sterile gauze pad and apply pressure.
5. Dispose of all used materials in appropriate container.




For Blood Collection Using A Syringe

Note: Should be performed wearing gloves.

1. Remove luer cap and attach syringe.
- Follow Blood Collection steps 2 – 5.



Note: Use the BD Vacutainer® Blood Transfer Device to transfer blood from the syringe to the evacuated tube. Dispose of all used materials in appropriate container.




For Short Term IV Administration (up to 2 hours)

Note: • Should be performed wearing gloves.

- When securing the device, be aware that taping over the button and applying pressure may cause the needle to retract.
- Prime set if necessary in accordance with your facility's procedure.

1. Remove needle sheath.
2. (a) Insert needle in vein by holding wings. (b) Observe for the presence of blood in the chamber.
3. Remove luer cap. Attach IV line or syringe to the female luer. Administer medication in accordance with your facility's procedure. Ensure that there is no air in the system.
 - The device should be changed in accordance with your facility's procedure and should not be used longer than 2 hours.
4. Retract the Needle: Depress the button. The needle will slide out of the venipuncture site and lock into place. Do not impede device retraction. Cover the puncture site with a sterile gauze pad and apply pressure.
5. Dispose of all used materials in appropriate container.



Priming Volume		
Gauge	Tubing (inch)	Volume (mL)
21	7	0.228
23	7	0.225
25	7	0.224
21	12	0.353
23	12	0.350
25	12	0.349

General Precautions

Handle all biologic samples and blood collection "sharps" (lancets, needles, luer adapters, and blood collection sets) in accordance with the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (e.g., through a puncture injury) since samples may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Utilize any safety engineered feature if the blood collection device provides one. Discard all blood collection "sharps" in biohazard containers approved for their disposal.

Using a syringe to draw blood may result in increased hemolysis.

GENERAL PRECAUTIONS: Read directions carefully before use.

Always follow Universal Precautions.

Activation of the device while the needle is still in the venipuncture site is recommended. Activation of the device after the needle is removed from the site should be performed away from self and others, since external blood droplets/IV fluid droplets may result.

Do not obstruct needle during activation. Visually confirm that the needle point is completely covered. After single-use, dispose of product according to your facility's regulations.

8364656
500002290
03/2013

Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417 USA
Becton, Dickinson and Company, Belliver Industrial Estate, Belliver Way,
Rohborough, Plymouth PL6 7BP United Kingdom
U.S. Patent Nos. D462,901, D465,571, D469,529 and 7,294,118. Made in USA
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2013 BD

STERILE R
Rx Only (USA)



Technical Questions?
www.bd.com/vacutainer/contact/
Call Toll-Free 1.800.631.0174 (USA Only)

Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417 USA
Becton, Dickinson and Company, Belliver Industrial Estate, Belliver Way,
Rohborough, Plymouth PL6 7BP United Kingdom
U.S. Patent Nos. D462,901, D465,571, D469,529 and 7,294,118. Made in USA
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2013 BD