**ConMed Linvatec 10k® / 24k® Pump Arthroscopy
Tubing Sets**

REF	Description
10k100	10k®/24k® Arthroscopy Inflow Tube Set
10k100D	10k®/24k® Arthroscopy Inflow Tube Set-Dyonics
10k150	10k® Arthroscopy Inflow/Outflow Tube Set

DESCRIPTION

These Tubing Sets are latex free and are provided sterile, single-use and are designed to be used with the ConMed Linvatec 10k and 24k Pumps.

INDICATIONS FOR USE

These Tubing Sets are intended for use, in conjunction with the ConMed Linvatec 10k and 24k Pumps, for pump-induced or gravity delivery of sterile fluid for irrigation and/or distention during arthroscopic procedures.

CONTRAINDICATIONS

- These Tubing Sets are not intended for hysteroscopic procedures.
- These Tubing Sets are not intended for pumping air or any other gas.


WARNINGS

- These Tubing Sets are single-use only. Do not clean, disinfect, resterilize or reuse. After use, dispose of properly.
- The ability to effectively clean and re-sterilize this single use device has not been established and subsequent re-use may adversely affect the performance, safety and/or sterility of the device.
- This device contains the phthalate plasticizer DEHP, which may leach into the patient during administration of irrigation. Animal studies have shown that high-level exposure to DEHP may affect male fertility and reproductive development; these findings have not been confirmed in clinical studies. In determining use of this device in a child, or a pregnant or nursing female, the clinical benefits should outweigh any risk to the infant.
- Avoid abrupt changes in joint position which may result in high intra-articular pressure spikes.
- Do not crimp, kink, puncture or roll over the tubing set with carts or other equipment. Erroneous reading may result, possibly causing injury and/or extravasation to the patient.
- Periodically examine the location of the inflow cannula to ensure that the cannula is within the joint capsule. Extravasation may occur as a result of improper cannula placement or excessive intra-articular joint pressure.

PRECAUTIONS

- Carefully read the ConMed Linvatec 10k and / or 24k Pump Instruction Manual for detailed instructions regarding the proper use and set-up of the pump console, pump tubing and pump accessories before using the ConMed Linvatec 10k or 24k Pump System.
- Maintain outflow when pump is in use.
- Frequently examine and palpate the patient's joint to check for extravasation and proper distension.

TUBING SET CONNECTION

- Install the appropriate Tubing Set using aseptic technique:
 - Circulator** - Open the tubing package. Present the package contents to the Scrub.
 - Scrub** - Remove the Tubing Set from its sterile package and pass the cassette and bag spike line out of the sterile field to the Circulator.
 - Circulator** - Position the tubing loop over the pump rotor with the ConMed Linvatec name on the cassette facing upward.
 - 10k: Clear plastic rotor
 - 24k: Central blue rotor
 - Circulator** - Pull the cassette down and stretch the tubing over the pump rotor. Position the cassette under the pump rotor and allow the cassette to snap into place.
 - Circulator** - Close the bag clamps.
 - Circulator** - Hang the fluid bag(s) on the infusion pole and spike the bags.
- Prime the Tubing Set:
 - Scrub** - Ensure the inflow clamp is open. Hold the patient end of the tubing over a basin.
 - Circulator** - Open the bag clamps.
 - Circulator** - Press the Run/Stop button on the pump or (Scrub Only) the remote control.
 - After all air is purged from the tubing, press the Run/Stop button again.
- Scrub** - Attach the tubing set to the inflow device.

Note: The pump display will indicate the SET pressure (mmHg) once the Tubing Set is connected.

10k150 Tubing Set Only

- Scrub** - Separate the outflow line from the inflow line, as needed, by pulling the two lines apart. Connect the outflow tubing line to the outflow device or shaver.
- Circulator** - Connect the suction / gravity tubing line to the appropriate hospital discharge mechanism.

PACKAGING and LABELING

- Tubing sets should only be used if the original packaging and labeling are intact.
- If packaging has been opened or altered, do not use. Sterile integrity may be compromised. Contact your local ConMed Linvatec representative, or in the U.S. contact the Customer Service Department.

CONTACT INFORMATION

For more information or a product demonstration, contact your local ConMed Linvatec sales representative, or call 1-800-237-0169 in the U.S.

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- Circulator** - Pull the cassette down and stretch the tubing over the pump rotor. Position the cassette under the pump rotor and allow the cassette to snap into place.

- Circulator** - Close the bag clamps.
- Circulator** - Hang the fluid bag(s) on the infusion pole and spike the bags.

- Prime the Tubing Set:
 - Scrub** - Ensure the inflow clamp is open. Hold the patient end of the tubing over a basin.
 - Circulator** - Open the bag clamps.
 - Circulator** - Press the Run/Stop button on the pump or (Scrub Only) the remote control.

- After all air is purged from the tubing, press the Run/Stop button again.

- Scrub** - Attach the tubing set to the inflow device.

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