



de Allgemeiner Benutzerhinweis

en General user information

fr Informations générales pour l'utilisateur

sv Allmän användarinformation

pt Informações gerais para utilizadores

nl Algemene gebruiksinstructies

da Generel brugerinformation

it Informazioni generali per l'utilizzatore

no Generell brukerinformatjon

es Información general para el usuario

en

General user information for single use products in sterile packaging

Please read this information carefully before use.

Warnings and cautions

- This product is to be used by qualified medical personnel only
- This product is for single use only. Do not re-sterilise and/or re-use
- Do not use if the package is open or damaged
- Re-tighten closure plugs before use
- Do not use if protection caps are displaced
- Stress cracks could appear when used for a longer time together with solutions based on alcohol and lipids
- For administration of blood and blood components, check suitability of infusion pump before application
- For disconnectable connections, check tightness before use
- If equipped with an orange protection cap (vented), replace it with a closure plug (non-vented)
- Observe the medication manufacturer's instructions for use regarding incompatibility and the present Summary of Product Characteristics
- When applied for infusion, use fluid filter <math>< 20 \mu\text{m}</math> (according to EN ISO 8536)
- Accepted good hygiene measures and working practices should be followed at all times
- The used product has to be disposed of correctly. Specially marked containers should be used necessary. Local regulations must be observed

Re-use warnings

The re-use of the product may cause unpredicted health risks! A used product is a potential carrier of pathogens, which means a health hazard to the patient and the spreading of infectious germs within the health care facility. Reprocessing of used products compromises safety. The manufacturer cannot be held responsible for the sterility assurance of re-used products within the health care facility.

P Pressure Compatible for use under pressure and gravity - according to EN ISO 1135 & EN ISO 8536	G Gravity Compatible for use under gravity only - according to EN ISO 1135 & EN ISO 8536
Non-pyrogenic Indicates a medical device that is non-pyrogenic. EN ISO 15223	XXX Contains or presence of Indicates the presence of identified product or substance. ISO 7000
Do not use if package is damaged Indicates a medical device that should not be used if the package has been damaged or opened. EN ISO 15223	XXX Does not contain or no presence of Indicates the no presence of identified product or substance. ISO 7000/EN 80416
XX ml Drops per millilitre Indicates the number of drops per millilitre. EN ISO 15223	PHT XXX Contains or presence of Phthalates Indicates the presence of identified phthalates. EN 15086
XX μm Liquid filter with pore size Indicates an infusion or transfusion system of the medical device that contains a filter of a particular nominal pore size. EN ISO 15223	One-way valve Indicates a medical device with a valve that allows flow in only one direction. EN ISO 15223
Do not re-use Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. EN ISO 15223	Pump, liquid pump ISO 7000
REF Catalogue number Indicates the manufacturer's catalogue number so that the medical device can be identified. EN ISO 15223	NON-STERILE Non-sterile Indicates a medical device that has not been subjected to a sterilization process. EN ISO 15223
LOT Batch code Indicates the manufacturer's batch code so that the batch or lot can be identified. EN ISO 15223	STERILE EO Sterilized using ethylene oxide Indicates a medical device that has been sterilized using ethylene oxide. EN ISO 15223
Date of manufacture Indicates the date when the medical device was manufactured. EN ISO 15223	STERILE R Sterilized using irradiation Indicates a medical device that has been sterilized using irradiation. EN ISO 15223
Manufacturer Indicates the medical device manufacturer, as defined in EU Directives 90/269/EEC, 93/42/EEC and 90/79/EEC. EN ISO 15223	Use-by date Indicates the date after which the medical device is not to be used. EN ISO 15223

Humidity limitation Indicates the range of humidity to which the medical device can be safely exposed. (e.g. between 30% and 90%) EN ISO 15223	Fragile, handle with care Indicates a medical device that can be broken or damaged if not handled carefully. EN ISO 15223
Temperature limit Indicates the temperature limits to which the medical device can be safely exposed. (e.g. between 15°C and 25°C) EN ISO 15223	Keep away from sunlight Indicates a medical device that needs protection from light sources. EN ISO 15223
Caution Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. EN ISO 15223	Keep dry Indicates a medical device that needs to be protected from moisture. EN ISO 15223
Consult instructions for use Indicates the need for the user to consult the instructions for use. EN ISO 15223	Language information (de = German, en = English, fr = French, sv = Swedish, pt = Portuguese, nl = Dutch, da = Danish, no = Norwegian, es = Spanish) ISO 639-1
Authorized representative in the European Community Indicates the authorized representative in the European Community. EN ISO 15223	I.V. For intravenous use only I.V. = latin: into "into" and versa "vein".
Filling/storage volume Tube volume is defined in millilitre per meter with an optional pressure definition. (e.g. 7 ml/1m (40°C) 9 ml/1m (40°C) 2 bar)) VOL	Not for intravenous use I.V. = latin: into "into" and versa "vein". IN
AirStop Indicates a medical device with a special AirStop filter membrane. AIR STOP	Do not use for more than 24 hours max. 24 h
Disconnectable connections disconnectable	Tube dimensions Indicates the tube dimension of the medical device. OD = Outer tube diameter. ID = Inner tube diameter (referring to the main line) OD (mm) ID