English

PRODUCT DESCRIPTION

Bionector/Bionector TKO is a catheter accessory, a needleless, closed IV system. It allows: sampling, injection and infusion. Bionector arterial, red coloured, is designed for arterial applications. Bionector/Bionector TKO can be used with all patient types. This device is not made with dry or natural rubber latex.

THREE TYPES OF PACKAGING

- rigid aseptic applicator (fig. 1)
- double-wrapped (rigid applicator and soft blister) (fig. 2)
- soft blister (fig. 3)

PRECAUTIONS OF USE

This device should be used by physicians and nurses who must know and apply all the standard precautions to prevent complications and problems related to the use of catheter accessories. Observe strict aseptic technique to prevent infections.

Sterility guaranteed unless unit pack has been damaged or opened.

To dispose of, place the contaminated items in the appropriate receptacle.

Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

Do not use any needle with Bionector/Bionector TKO (fig.4).

Do not put a cap, plug or obturator of any sort on Bionector/Bionector TKO other than a disinfecting cap (fig.5).

Do not use clamps or forceps to connect/tighten or disconnect the male Luer-lock administration set to Bionector/Bionector TKO, as this can lead to overtightening and damage to Bionector/Bionector TKO.

Never connect Bionector/Bionector TKO onto another Bionector/Bionector TKO. When Bionector/Bionector TKO is on the port of a ramp or a stopcock, it is recommended that the port should be in a closed position when not in use.

When Bionector/Bionector TKO is present on multi lumen extension tubes, it is recommended that the lumen should be clamped when not in use.

Bionector/Bionector TKO should be used with devices with a male Luer that conforms to ISO standards related to Luer connectors, and fitting whose internal diameter is equal or greater than 1.40 mm.

DISINFECTING BIONECTOR/BIONECTOR TKO

1. Disinfect Bionector/Bionector TKO before use and after use unless a disinfecting cap is in place (fig.6).
2. Clean in accordance with your hospital/department protocol.
3. Bionector/Bionector TKO can be cleaned with most disinfecting agents. However antiseptics containing both alcohol and hydrogen peroxide can damage Bionector/Bionector TKO and must be avoided.
4. Allow a minimum of 30 seconds for the disinfection agent to dry before use.
5. Do not enclose or wrap Bionector/Bionector TKO in materials or containers impregnated with disinfectants.

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








English

CONNECTION TO BIONECTOR/BIONECTOR TKO

You will find that making a Luer-lock connection on Bionector/Bionector TKO simply needs a little more pressure than traditional Luer fittings (fig.7).

When using a Luer Slip fitting, e.g. a non-locking syringe, you will find that a quarter turn of the syringe as it is pushed firmly in will help provide a secure fit (fig.8).

When using Bionector wrapped in a rigid protective non-touch applicator, connect the female Luer to the Bionector in its applicator, flex it up and down until you hear a crack to release Bionector from its applicator (fig. 9).

F.A.Q.

Q. How often do I need to replace Bionector/Bionector TKO?

A. Every 7 days or 360 accesses.

Q. What is the flow rate of Bionector/Bionector TKO?

A. The flow rate of Bionector is 105mL/min. The flow rate of Bionector TKO is 55mL/min [gravity feed ISO 10555].

Q. Can I give blood or blood products via Bionector/Bionector TKO?

A. Yes.

Q. Can I infuse lipid or lipid products through Bionector/Bionector TKO?

A. Yes, Bionector/Bionector TKO is made of a lipid resistant material.

Q. Is Bionector/Bionector TKO compatible to use with MRI scanners?

A. Yes, Bionector/Bionector TKO is conditional for use with MRI scanners. However MR images may be compromised if the area of interest is within 10 cm proximity to the device; in this case, artefacts are likely to occur. To avoid this, simply move Bionector/Bionector TKO approximately 10 cm from the area being-scanned.

Q. Do I need to prime Bionector/Bionector TKO?

A. In the case of the range of short extensions with Bionector/Bionector TKO, yes. In the case of the individual Bionector, the priming volume is 0.03 mL and thus negligible from a clinical viewpoint. However there may be specific occasions when you may wish to prime Bionector. For Bionector TKO, the priming volume is 0.1mL, therefore we recommend to prime it.

Q. Is Bionector/Bionector TKO CT-Rated?

A. Yes, for the single-bung Bionector/Bionector TKO. Maximum pressure rated to 350 psi (24 bars).

For delivery of the contrast media (Visipaque 320 mg Iodine/mL), the flow rate was measured as below:

Infusion pressure	Flow rate mL/sec
35 psi (2.4 bars)	> 5 mL/s
145 psi (10bars)	> 10 mL/s

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