



Edge™

Electrodes

REF E1455
REF E1455-4
REF E1455-6
REF E1465
REF E1465-4
REF E1465-6

en Instructions for Use
fr Mode d'emploi
pt Instruções de uso
da Brugervejledning
de Gebrauchsanweisung
es Instrucciones de uso
fi Käyttöohjeet
it Istruzioni per l'uso
nl Gebruiksaanwijsties
no Bruksanvisning
sv Bruksanvisning
ru Инструкции по применению
zh 使用说明

REF E1455

Insulated Blade Electrode
2.75" (7.0 cm)

REF E1455-4

Insulated Blade Electrode
4" (10.16 cm)

REF E1455-6

Insulated Blade Electrode
6.5" (16.51 cm)

REF E1465

Insulated Needle Electrode
2.8" (7.2 cm)

REF E1465-4

Insulated Needle Electrode
4" (10.16 cm)

REF E1465-6

Insulated Needle Electrode
6.5" (16.51 cm)



Not made with natural rubber latex



Do not use if package is opened or damaged

For use with a max peak voltage of 5600 V.

Indications for Use

The Edge Coated Electrodes are intended as an alternative to uncoated stainless steel electrodes for use in conventional monopolar electrosurgical accessories. The coated electrodes are intended for use in situations where monopolar electrosurgical cutting and coagulation are normally used.

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in bio-incompatibility, infection, or product failure risks to the patient.

Warning

Danger: Explosion Hazard Do not use electrosurgery in the presence of flammable anesthetics.

Fire Hazard Both oxygen (O₂) and nitrous oxide (N₂O) support combustion. Watch for enriched O₂ and N₂O atmospheres near the surgical site, especially during head and neck surgery. Enriched O₂ atmospheres may result in fires and burns to patients or surgical personnel.

Warning

Fire/Explosion Hazard The following substances contribute to increased fire and explosion hazards in the operating room:

- Oxygen-enriched environments
- Oxidizing agents, such as N₂O atmospheres
 - Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
 - Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
 - If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.
 - If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.
- Alcohol-based skin prep/ping agents and tinctures
 - Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
- Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.

Warning

Fire Hazard The sparking and heating associated with electrosurgery can provide an ignition source.

Observe fire precautions at all times:

- When using electrosurgery in the same room with gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes or near the surgical site.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen-enriched environments, such as in throat or mouth procedures. Eschar plus high oxygen may create embers. Keep the electrode clean and free of all debris.
- Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

The electrode must fit completely and securely into the pencil. An incorrectly seated electrode may result in burns to the patient or surgical personnel.

4

Warning

Fire Hazard Always place the active electrode in a clean, dry, insulated safety holster when not in use.

- Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.
- Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials, such as gauze or surgical drapes. Place longer electrodes, such as extended electrodes, away from the patient and drapes.

Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

Warning

Some surgeons may elect to "buzz the hemostat" during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon's hands may result. To minimize the risk, take these precautions:

- Do not "buzz the hemostat" with a needle electrode.
- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate cut rather than coag. Cut has a lower voltage than coag.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.
- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes the current concentration at the fingertips.

Warning

"Buzz the hemostat" below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon's hands.

- When using a coated or nonstick blade electrode, place the **edge** of the electrode against the hemostat or other metal instrument.


Precaution

Do not exceed maximum power limits as stated in instructions for use. Exceeding these power settings may result in patient injury or product damage.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.


5

2. Grasp the insulating sleeve on the electrode and insert the electrode into the pencil.




3. Ensure the insulating sleeve fits securely inside the nose of the pencil so that the nose overlaps the insulating sleeve by at least 1/8" (0.3 cm).

4. Hex electrodes have a depth indicator. The line on the depth indicator should be flush with the tip of the handswitching pencil.



5. A tip protector covers the coated end of some electrodes. If a tip protector is present, remove it before use.



Maximum Power

Maximum Power for Blades

The maximum power limits for E1455 series blades are as follows:

Mode	Power Setting
Coag	35 watts
Pure Cut or Blend	50 watts

Maximum Power for Needles

The maximum power limits for E1465 series needles are as follows:

Mode	Power Setting
Coag	25 watts
Pure Cut or Blend	30 watts

Part No. PT00130893

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.™™™ brands are trademarks of their respective owner.

STERILE EO



Single use

**RX
ONLY**



Consult
instructions
for use



Do not re-sterilize





MEZO



2797

© 2012 Covidien.

 Manufactured for Covidien llc,
15 Hampshire Street, Mansfield, MA 02048 USA.
 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made in Mexico. Printed in Mexico.

www.covidien.com

REV 08/2020