

CODAN



- de** Allgemeiner Benutzerhinweis
- en** General user information
- fr** Informations générales pour l'utilisateur
- sv** Allmän användarinformation
- pt** Informações gerais para utilizadores
- nl** Algemene gebruiksinstructies
- da** Generel brugerinformation
- it** Informazioni generali per l'utilizzatore
- no** Generell brukerinformasjon
- es** Información general para el usuario

en

General user information

for single use products in sterile packaging

Please read this information carefully before use.

















Warnings and cautions

- This product is to be used by qualified medical personnel only
- This product is for single use only. Do not re-sterilise and/or re-use
- Do not use if the package is open or damaged
- Re-tighten closure plugs before use
- Do not use if protection caps are displaced
- Stress cracks could appear when used for a longer time together with solutions based on alcohol and lipids
- For administration of blood and blood components, check suitability of infusion pump before application
- For disconnectable connections, check tightness before use
- If equipped with an orange protection cap (vented), replace it with a closure plug (non-vented)
- Observe the medication manufacturer's instructions for use regarding incompatibility and the present Summary of Product Characteristics
- When applied for infusion, use fluid filter <math>< 20 \mu\text{m}</math> (according to EN ISO 8536)
- Accepted good hygiene measures and working practices should be followed at all times
- The used product has to be disposed of correctly. Specially marked containers should be used if necessary. Local regulations must be observed

Re-use warnings

The re-use of the product may cause unpredicted health risks! A used product is a potential carrier of pathogens, which means a health hazard to the patient and the spreading of infectious germs within the health care facility. Reprocessing of used products compromises safety. The manufacturer cannot be held responsible for the sterility assurance of re-used products within the health care facility.

 Pressure	Compatible for use under pressure and gravity - according to EN ISO 1135 & EN ISO 8536	 Gravity	Compatible for use under gravity only - according to EN ISO 1135 & EN ISO 8536
 EN ISO 15223	Non-pyrogenic Indicates a medical device that is non-pyrogenic.	 ISO 7000	Contains or presence of Indicates the presence of identified product or substance.
 EN ISO 15223	Do not use if package is damaged Indicates a medical device that should not be used if the package has been damaged or opened.	 ISO 7000/EN 80416	Does not contain or no presence of Indicates the no presence of identified product or substance.
 EN ISO 15223	Drops per millilitre Indicates the number of drops per millilitre.	 XXX EN 15986	Contains or presence of Phthalates Indicates the presence of identified phthalate.
 EN ISO 15223	Liquid filter with pore size Indicates an infusion or transfusion system of the medical device that contains a filter of a particular nominal pore size.	 EN ISO 15223	One-way valve Indicates a medical device with a valve that allows flow in only one direction.
 EN ISO 15223	Do not re-use Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	 ISO 7000	Pump, liquid pump
 EN ISO 15223	Catalogue number Indicates the manufacturer's catalogue number so that the medical device can be identified.	 EN ISO 15223	Non-sterile Indicates a medical device that has not been subjected to a sterilization process.
 EN ISO 15223	Batch code Indicates the manufacturer's batch code so that the batch or lot can be identified.	 EN ISO 15223	Sterilized using ethylene oxide Indicates a medical device that has been sterilized using ethylene oxide.
 EN ISO 15223	Date of manufacture Indicates the date when the medical device was manufactured.	 EN ISO 15223	Sterilized using irradiation Indicates a medical device that has been sterilized using irradiation.
 EN ISO 15223	Manufacturer Indicates the medical device manufacturer, as defined in EU Directives 90/386/EEC, 93/42/EEC and 98/79/EC.	 EN ISO 15223	Use-by date Indicates the date after which the medical device is not to be used.

 EN ISO 15223	Humidity limitation Indicates the range of humidity to which the medical device can be safely exposed. (e.g. between 30% and 60%)	 EN ISO 15223	Fragile, handle with care Indicates a medical device that can be broken or damaged if not handled carefully.
 EN ISO 15223	Temperature limit Indicates the temperature limits to which the medical device can be safely exposed. (e.g. between 15°C and 25°C)	 EN ISO 15223	Keep away from sunlight Indicates a medical device that needs protection from light sources.
 EN ISO 15223	Caution Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	 EN ISO 15223	Keep dry Indicates a medical device that needs to be protected from moisture.
 EN ISO 15223	Consult instructions for use Indicates the need for the user to consult the instructions for use.	 ISO 639-1	Language information (de = German, en = English, fr = French, sv = Swedish, pt = Portuguese, nl = Dutch, da = Danish, it = Italian, no = Norwegian, es = Spanish)
 EN ISO 15223	Authorized representative in the European Community Indicates the authorized representative in the European Community.		For intravenous use only I.V. = latin: intra "into" and vena "vein".
 EN ISO 15223	Filling/storage volume Tube volume is defined in millilitre per meter with an optional pressure definition. (e.g. 7 ml/1 m (40°C) 9 ml/1 m (40°C/2 bar))		Not for intravenous use I.V. = latin: intra "into" and vena "vein".
	AirStop Indicates a medical device with a special AirStop filter membrane.		Do not use for more than 24 hours
 disconnectable	Disconnectable connections		Tube dimensions Indicates the tube dimension of the medical device. OD = Outer tube diameter, ID = Inner tube diameter (referring to the main line)