smiths medical

portex®

Spinal Needle Systems

en Instructions for use

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sv Bruksanvisning

da Brugsanvisning

no Bruksanvisning

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www.smiths-medical.com

en **CE** **2797**

en **!** **Caution** **REF** Catalog Number **!** Not made with natural rubber latex **!** Do Not Reuse **!** Date of Manufacture **!** Do not use if package is damaged **!** Keep dry **!** Keep away from sunlight **LOT** Batch Code **!** Use By **!** Manufacturer **!** Quantity **STERILE EO** Sterilized using ethylene oxide **EC REP** Authorized Representative in the European Community

en **!** **!** **!** **!** **STERILE EO** **REF** **LOT** **!** **!** **!** **!** **EC REP** **!** **!**

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en *English*

Portex® Spinal Needle Systems

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Instructions for use

- Check that the spinal needle will pass freely through the introducer needle.
- Perform aseptic preparation of the skin.
- Infiltrate local anaesthetic into the skin, subcutaneous and deeper tissues.
- Insert the introducer needle, avoiding excessive insertion which might pierce the dura. Check that the stylet of the spinal needle is correctly fitted. Insert the spinal needle through the introducer needle, and advance it through the resistances of first the ligamentum flavum, and then the dura mater. After piercing the dura, the needle should be rotated so that the lateral eye (pencil point) is less likely to be occluded by a nerve root.
- Withdraw the stylet from the spinal needle and wait for the cerebrospinal fluid (CSF) to appear in the needle hub. Attach the syringe containing local anaesthetic and aspirate gently to confirm free flow of CSF. If CSF does not appear in the needle hub, replace the stylet and rotate the needle 90° or advance or withdraw the needle until CSF appears and confirm free flow of CSF in the syringe.
- Inject the local anaesthetic, withdraw the spinal and introducer needles, and cover the puncture site with a suitable dressing.

Precautions

- Ensure that the bevel on lancet point needles is facing laterally during puncture of the dura.
- Always advance or re-position the spinal needle with the stylet in place.
- If there is any resistance to injecting the local anaesthetic through the spinal needle, or pain on injection, stop immediately and proceed in accordance with currently accepted medical techniques.

Warnings

- Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.**
- Never use excessive force to advance the spinal needle when bone is felt. Continued advancement may bend the tip of the spinal needle. If undue resistance is felt, carefully remove the spinal needle and discard.**
- If resistance is encountered when removing or re-positioning the spinal needle, the spinal needle and introducer needle must be withdrawn together and the procedure repeated using new needles.**
- Always use spinal needles with an introducer or Tuohy needle to reduce the risk of the spinal needle bending. The tip of a bent spinal needle may break inside the patient when attempt is made to withdraw it. A bent spinal needle, or the surgical intervention required to remove a broken needle tip, will cause unnecessary trauma to the patient.**

Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2017, 2019 Smiths Medical. All rights reserved.

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