ETHICON™

BONE WAX

- STERILE -

DESCRIPTION

ETHICON™ Bone Wax is a sterile mixture of beeswax, paraffin, and isopropyl palmitate, a wax-softening agent. It is opaque and has a waxy odor.

INDICATIONS

Bone Wax is intended for use for the control of bleeding from the divided, drilled or chipped edges of bone by physically plugging the osseous canals which contain the bleeding capillaries.

ACTIONS/PERFORMANCE

Bone Wax achieves local hemostasis of bone by acting as a mechanical (tamponade) barrier. It does not act biochemically and is nonabsorbable.

CONTRAINDICATIONS

Bone Wax should not be used where rapid osseous regeneration and fusion are desired.

WARNINGS/PRECAUTIONS/INTERACTIONS

Bone Wax should not be subjected to excessive heat.

Bone Wax should be used sparingly and any excess removed from the operative site.

Bone Wax may inhibit osteogenesis and may act as a physical barrier to the reparative process.

The package containing Bone Wax should be opened just prior to use to minimize the potential of contamination and excessive drying.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

ADVERSE REACTIONS

Adverse reactions associated with the use of this device include mild inflammatory reactions in tissues adjacent to the site of implantation and inhibition of osteogenesis. Like all foreign bodies, Bone Wax may potentially cause an existing infection.

DOSAGE ADMINISTRATION

Bone Wax should be used immediately after removal from the package. Using aseptic technique, Bone Wax should be warmed to desired consistency by manipulation with the fingers. The softened Bone Wax is then applied to the bone edges as indicated by surgical circumstances and the preference of the surgeon.

HOW SUPPLIED

Bone Wax is available sterile in individual foil envelopes, each containing 2.5 grams, and packaged in an individually sealed overwrap packet.

STERILITY

Bone Wax is sterilized by irradiation. Do not resterilize. Do not use if package is opened or damaged. Discard opened, unused Bone Wax.

STORAGE

No special storage conditions required. Do not use after expiry date.

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