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**IMPORTANT INFORMATION**  
Please Read Before Use

**CODMAN® Surgical Patties**  
**CODMAN Surgical Strips**

**FRANÇAIS**

**STERILE** **EO**

**Description**  
CODMAN® Surgical Patties and CODMAN Surgical Strips are manufactured of COTTONOID® Material with x-ray detectable markers. All patties have a suture string attached for ease in performing post-surgical count verification.

**Indications**  
The surgical patties and surgical strips are indicated for the use in protection of tissue, including brain and other tissues of the central nervous system, during surgery.

**Warnings**  
Do not leave patties or strips in situ. Failure to remove them may result in a foreign body reaction.

Avoid cutting the patties or strips because fragments without x-ray detectable material may enter the surgical site. Fragments left in the surgical site may result in an unintended adverse reaction, including a foreign body reaction.

**Precautions**  
The string attached to the patties is for identification purposes only. Do not use the string to pull the patty from the wound.

This device is considered radiopaque in accordance with ASTM Test Method F640 (Standard Test Methods for Radiopacity of Plastics for Medical Use). The device may not be visible under all radiographic conditions. The size and position of the marker may impact radiopacity. In addition, adjacent structures may inhibit visibility.

It is recommended that at least three views, using the optimal parameters for the imaging (x-ray) equipment, at a variety of angles (e.g., 45 degrees, 22.5 degrees, and 0 degree angles) for anterior and posterior, or the appropriate plane, be taken and examined for a missing device. If there are

concerns regarding visualization, consult with your local imaging expert to establish the optimal radiographic parameters (e.g., kVp, mAs) for visualization with the imaging equipment.

The locator strings attached to the patties will not be visible by x-ray as they are not radiopaque. The locator strings are for identification purposes only; avoid using the strings to remove the patties from the surgical site to prevent detachment.

Count all devices before and after the procedure prior to surgical closure. In the event a device cannot be located, an x-ray can be used to locate the device. Only the radiopaque markers are visible on imaging. The size and position of the radiopaque markers may impact their visibility.

**Sterility**  
②  
This product is for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Integra Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. These devices are intended to come into contact with the central nervous system and the ability does not currently exist to destroy possible contaminates such as Creutzfeldt-Jakob Disease. Resuse can also compromise device performance and any usage beyond the design intent of this single-use device can result in unpredictable use hazards or loss of functionality.

Integra LifeSciences will not be responsible for any product that is resterilized, nor will we accept for credit or exchange any product that has been opened but not used.

As long as the individual package of patties or strips is not opened or damaged, the product is sterile.

**Instructions for Use**  
Verify count of patties or strips when package is opened. Moisten the patties or strips with sterile water or physiologic saline before applying to tissue. Replace patties on the count card as they are removed to facilitate count verification.

**CAUTION:** Failure to perform counts of patties or strips may result in additional procedures, prolonging surgery.

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**PRODUCT INFORMATION DISCLOSURE**  
INTEGRA LIFESCIENCES CORPORATION ("INTEGRA") HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. INTEGRA WARRANTS THAT THESE PRODUCTS SHALL CONFORM TO THE PRODUCT LIMITED WARRANTY AS PROVIDED IN THE PRODUCT LABELING OR APPLICABLE PRODUCT CATALOG. THIS WARRANTY IS EXCLUSIVE, AND INTEGRA DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THESE PRODUCTS. INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

Codman and Cottonoid are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.

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