
# Spinocan® – Spinal Needle for Single Shot Technique

---

---

**B | BRAUN**

LLDorder 5899 - G 220799

<b>GB</b>	Spinocan® – Spinal Needle for Single Shot Technique	<b>6</b>
<b>ES</b>	Spinocan® – Aguja espinal para técnica de inyección única	<b>8</b>
<b>FR</b>	Spinocan® – Aiguille spinale pour technique single shot	<b>10</b>
<b>ID</b>	Spinocan® – Jarum Spinal untuk Teknik Single Shot	<b>12</b>
<b>KR</b>	Spinocan® - 단발 기법용 척추 바늘	<b>14</b>

			<b>LOT</b>	<b>STERILE</b>
<b>GB</b>	Do not re-use	Consult instruction for use	Batch number	Sterilized using ethylene oxide
<b>ES</b>	No reutilizable	Consulte las instrucciones de uso	Código de lote	Esterilizado con óxido de etileno
<b>FR</b>	Ne pas réutiliser	Consulter les instructions d'utilisation	Numéro de lot	Stérilisation à l'oxyde d'éthylène
<b>ID</b>	Jangan gunakan kembali	Baca petunjuk penggunaan	Nomor bets	Disterilkan menggunakan etilena oksida
<b>KR</b>	재사용금지	사용 지침서 확인	로트 번호	에틸렌 옥사이드 멸균

				
Use-by date	Manufacturer	Date of manufacture	Do not use if package is damaged	Caution
Fecha de caducidad	Fabricante	Fecha de fabricación	No utilizar si el envase está dañado	Atención
Date limite d'utilisation	Fabricant	Date de fabrication	Ne pas utiliser si l'emballage est abîmé	Attention
Gunakan sebelum	Produsen	Tanggal produksi	Jangan gunakan jika kemasan rusak	Perhatian
사용기한	제조업체	제조일	포장이 손상된 경우 사용하지 마십시오	주의

## **GB** Instructions for use

### Materials used

Stainless steel, nickel plated brass, polycarbonate, polypropylene, epoxid, resin

### Indications

- Lumbar puncture (diagnosis)
- Spinal anesthesia (Injection of local anesthetics into the subarachnoid space for diagnosis, operative procedures and pain management).

### Contraindications

The recognized contraindications for spinal anesthesia and lumbar puncture must be observed.

Be aware of:

- blood coagulation disorders
  - skin infection at or near the puncture / injection site
  - sepsis
  - History of hypersensitivity to local the anaesthetic or any of the materials employed is known
  - Patient refusal
  - Severe untreated hypovolaemia
  - Increased intracranial pressure
- Spinal anesthesia should not be performed in patients with the following conditions:
- severe decompensated hypovolemia
  - shock
  - acute cerebral or spinal cord disease
  - increased intracranial pressure
  - existing infection at the injection site

- endogenous or iatrogenic coagulopathy
- Anatomical deformities of the patient's back
- Amniocentesis

### Risks

The known risks of lumbar puncture and spinal anesthesia are especially post dural puncture headache, hypotension and bradycardia, cardiac arrest, apnoea or acute toxicity to the local anesthetic. Neurological disorders such as a consequence of epidural haematoma or abscess formation, anterior spinal artery syndrome, cauda equina syndrome, urinary retention, herniation, pain, infections, mispositioning, hypothermia, nausea and vomiting, total spinal anaesthesia, low cerebrospinal fluid pressure syndrome etc. are very rare but still recognized complications. Lumbar puncture and spinal anesthesia should only be performed with the appropriate equipment and trained staff to manage such complications.

Very rare complications:

- hearing loss
- epidermoid tumor
- retroperitoneal abscess

### Warning

**Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may**

**lead to injury, illness or death of the patient.**

**Do not use if packaging is damaged. Do not re-sterilize.**

**The needle might bend and eventually break when applying too strong forces during advancement. The tip of the needle becomes blunt or can even split after bone contact. Needles with damaged tip increase the risk of post dural puncture headache due to an extended perforation of the dura mater.**

### Duration of use

Spinal needles are used for a few minutes for diagnostic puncture (lumbar puncture) or for the injection of drugs into the subarachnoid space (spinal anesthesia) and are withdrawn from the patient immediately after completion of the procedure.

### Instructions for use

Lumbar puncture/spinal anesthesia

**Aseptic conditions are required.**

1. Orient the patient according to his/her condition and/or the planned puncture/injection site. Choose between lateral, sitting or prone position. Locate the midline and follow it during injection. Try to minimize lumbar lordosis.
2. Select a suitable spinal needle and use an Introducer Needle for spinal needles less than 25G diameter. Check emer-

- gency equipment (intubation kit, drugs) and ensure continuous monitoring of the patient.
3. Disinfect puncture site.
  4. Inject 1–1.5 cc of local anesthetic into the skin at the selected puncture site of the spinal needle to anesthetize the tissue and the supraspinous and intraspinal ligaments.
  5. The guiding needle is introduced into the puncture site until reaching the interspinous ligament. The guiding needle is fixed with thumb and index finger of the left hand (for right-hand user).
  6. The spinal needle including the stylet is advanced through the guiding needle. The typical "dural click" occurs when the ligamentum flavum has been passed and the subarachnoid space has been reached.
- Warning**  
If any resistance is felt during advancement of the spinal needle carefully correct the orientation of the needle but never apply strong forces in order to overcome the obstacles.
7. Remove the stylet from the spinal needle. Carefully check the occurrence of cerebrospinal fluid which can slightly be contaminated by blood but must never be pure blood.
    - Lumbar puncture procedure – (take liquor for diagnosis, withdraw spinal needle,

- cover the puncture site with sterile plaster).
8. – Spinal anesthesia – Before and shortly after injection of the local anesthetic withdraw liquid from the subarachnoid space. Inject the drug as described in the manufacturer's instructions.
- Note**
- The local anesthetic must never be injected without the evidence of cerebrospinal fluid.
  - If pure blood is flowing instead of cerebrospinal fluid withdraw the needle and re-try at a different intervertebral space.
  - If no cerebrospinal fluid is flowing back rotate the needle to all four quadrants and carefully aspirate until cerebrospinal fluid is visible. If this procedure does not show any backflow of cerebrospinal fluid repeat the spinal puncture using a different needle direction.
  - If paresthesia occurs after reaching the subarachnoid space the needle must be slightly withdrawn.
  - If paresthesia occurs during injection the needle must be repositioned before continuing the injection.
9. After complete injection of the drug withdraw the needle and cover the puncture site with

- sterile plaster.
10. Spreading of the anesthetic drug is done by turning the patient into a favorable position depending on the type of drug used (hyperbaric, isobaric or hypobaric technique).
  11. The level of local anesthesia is to be continued and carefully monitored by e.g. touching the skin with a cold Ecoflac bottle. After reaching complete block of the motoric and/or sensoric and/or sympathetic nerves in the desired dermatom, surgery can start.

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B. Braun Melsungen AG