en **Instructions for use - Infusion set**

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
STERILE EO MD

Infusion set

suitable for use with compatible pumps as well as gravity feed

Intended purpose

Transfer of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions by intravenous infusion.

Indications

Diseases and conditions requiring intravenous infusion of fluids, pharmaceuticals (excluding blood and blood components) and parenteral nutrition solutions.

Contraindications

The product is not designed or approved for:

- Home use
- Use by unqualified, untrained staff

Self-application is not permitted and is not considered as intended use.

The product is not intended for the transfer of:

- Blood and blood components
- Enteral nutrition solutions

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions with regard to incompatibilities (interaction material/ pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed. Special attention should be paid to vulnerable patient groups such as children and adolescents as well as pregnant or lactating women.

Patient target group

Patients of at least 2 years of age and/or a weight of 12 kg.

The product can be used for any adult patient and for permitted paediatric patient populations with consideration given to the adequacy of vascular anatomy, appropriateness of the solution being infused, and the duration of therapy.

User

Qualified medical and pharmaceutical personnel only.

For use in hospitals, specialised practices and pharmacies with conventional environmental and hygiene conditions.

Clinical benefits

The transfer of fluids, pharmaceuticals and parenteral nutrition solutions by intravenous infusion can bring about therapeutic benefits that improve the health of the patient.

Infusions of fluids, pharmaceuticals and parenteral nutrition solutions carry implicit risks and should therefore only be carried out if they improve the patient's health.

General safety instructions

- The product is a medical device. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product has to be completely disposed of and replaced by a new one.
- The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage, detached protective caps) or performance changes, the product shall not be used on the patient, and has to be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • Prior to the transfer of fluids, pharmaceuticals and parenteral nutrition solutions, the qualified medical personnel shall conduct a patient-specific benefit-risk assessment.
- Furthermore, the international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account. • Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • In case of consecutive administration of different drugs with the same infusion set (cf. exemplary product illustrations C and D), rinsing with a neutral compatible solution is recommended after each medication to avoid potential incompatibility of drug mixtures, maintain the integrity of the

line and prevent microbiologic growth. • The function of the product shall be checked by qualified medical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods have to be observed at all times. • Use of protective medical gloves is recommended to prevent any infections of users.

Notes on safe use of product

- The product is designed for use with compatible pressure infusion pumps with generated pressures up to 200 kPa (2 bar). The product can also be used under gravity and with a pressure cuff with generated pressures up to 50 kPa (0.5 bar). Restrictions in the type of application are possible due to the product structure (e.g. products with integrated pressure-activated valve cannot be used by gravity feed).
- Notes on use with pumps: the user shall ensure that the pump is approved for the administration of fluids, pharmaceuticals and parenteral nutrition solutions and calibrated with the appropriate infusion device before commissioning. Always observe the instructions for use from the pump manufacturer. Only open the roller clamp during pump operation!
- Before using the product, compare the permissible fluid intake of the patient with the suitability of the product. • Prior to use, the product shall be checked for damage, tightness and any misconnections.
- Prior to use, existing closure plugs and disconnectable connections have to be checked for tightness and retightened if necessary. • If equipped with an orange protective cap (vented), replace it with a closure plug (non-vented) before using the medical device. • A fluid filter < 20 µm (according to EN ISO 8536-8) is required for infusion.
- Priming of products equipped with FlowStop protection cap: prime the product only by gravity feed with a maximum pressure of 0.15 bar (respecting a maximum water column of 150 cm). Do not exert manual pressure on the infusion container. If higher pressure is applied on the FlowStop protection cap, its function cannot be assured (i.e. protection against leakages). In this case, the FlowStop protection cap has to be replaced by a sterile closure plug immediately after priming. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70 % or isopropyl alcohol 70 % in combination with chlorhexidine gluconate 2 % can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • Long-term administration of alcohol- and lipid-based solutions can cause stress cracks. • Excessive forces on Luer-Lock connections can lead to stress cracks. • Avoid excessive injection pressure during injections of solutions into the infusion device/ line, especially when using small syringes. Handling errors can result in leakages.

Use of product

Please also observe the exemplary product illustration(s) as well as the graphic brief instructions under notes on product use and all further information/images for product use on and in the product packaging.

General use (cf. exemplary product illustrations A / B / C / D):

1. Check if the sterile packaging is damaged or has been unintentionally opened.
2. Remove the product from the sterile packaging and check for possible damage.
3. In case of rigid/semi-rigid infusion containers/bottles, an infusion device with a vented closure-piercing device (1) (A / D) is required. In case of infusion bags, infusion devices with non-vented (1) (B / C) as well as vented closure-piercing devices (1) (A / D) can be used. Close the ventilation cap of vented closure-piercing devices (1) (A / D) first.
4. Close the roller clamp (6).
5. Remove the protective cap and insert the vented closure-piercing (1) (A / D) straight into the uprightly positioned rigid/semi-rigid infusion container/bottle or insert the non-vented closure-piercing device (1) (B / C) into the port of infusion bag with a rotating movement. Please ensure that the wall of the bag remains intact during piercing.
6. Fill the drip chamber (5) up to the level ring (4).
7. In the case of rigid/semi-rigid infusion containers/bottles, open the ventilation cap of the vented closure-piercing device (1) (A / D). For infusion bags, the ventilation cap remains closed.
8. Open the roller clamp (6) and prime the remaining infusion device up to the patient connection port (9) so that the product is applied without any air inclusions. Close the roller clamp (6).

9. In case of gravity application: proceed according to point 10.
When using with pump: follow the pump manufacturer's instructions for use. Open the door of the infusion pump. Insert the tubing into the infusion pump. Roller clamp (6) has to be positioned on the patient connection port (9) side. Close the door of the infusion pump.
10. Connect the patient connection port (9) to the patient access.
11. For starting the infusion, open the roller clamp (6) and adjust the required medication and patient parameters.
12. Check the infusion at regular intervals.

Use of product configurations with connection port above/integrated into the drip chamber (cf. exemplary product illustration C, e.g. Drip SWAN®) for secondary line(s)/single-use syringe(s):

A secondary line for the administration of an infusion container or a single-use syringe that has been pre-filled for infusion can be connected to the disinfected connection port (3) above/integrated into the drip chamber (5) of the infusion device. Successive administration of infusion solutions through several secondary lines or single-use syringes is possible. The respective unused line shall be closed with the corresponding closure clamp. After each medication administration, flushing with a neutral compatible solution is recommended.

Use of product configurations with multiple closure-piercing devices (cf. exemplary product illustration D):

The various closure-piercing devices (1) can be used for successive administration of infusion solutions. The respective unused line(s) shall be closed with the corresponding closure clamp(s) (2). After each medication administration, flushing with a neutral compatible solution is recommended.

Regarding products equipped with Wing-Valve drip chamber (5) (D), the drip chamber can be filled or regulated in the liquid level by pressing the Wing-Valve wings together.

Product compatibility

• The closure-piercing device complies with EN ISO 8536-8 and can be connected with all compatible infusion containers. • All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/

EN ISO 80369-7. • Compatible with volumetric infusion pumps with a pressure alarm limit not exceeding 2.0 bar.

Residual risks

• Miscalculations of the Luer connectors with medical devices of other application areas (e.g. enteral nutrition) are possible. • In the event of handling errors or, in rare cases, material- and/or production-related limitation of product functionality, leakages can occur, potentially resulting in air embolism and/or under-/overdosing of the infusion. • The aforementioned residual risks can lead to a significant deterioration of the patient's health status and possibly even to death depending on the patient's health condition.

Period of use

The international and national guidelines for the transfer of fluids, pharmaceuticals and parenteral nutrition solutions shall be taken into account.

A maximum period of use of 96 hours should not be exceeded.

A maximum period of use of 96 hours is recommended for non-lipid-containing solutions and 24 hours for lipid-containing solutions.

Taking into account the individual benefit-risk ratio of each patient and the Summary of Product Characteristics of the pharmaceutical manufacturer, a longer period of use can be considered by the physician.

Product disposal

• The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. • The closure-piercing device can be inserted into the notch on the roller clamp to avoid any injuries. • Protective medical gloves have to be worn during disposal.

Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay.

Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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