## **en** Instructions for use

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<https://www.codancompanies.com/download-centre/>

## Transfer set for blood and blood components

### Withdrawal spike for blood bags TAKE SET B93/TAKE SET B88

#### Intended purpose

Extracorporeal withdrawal, preparation (mixing, reconstitution) or transfer of blood and blood components and, if necessary, electrolyte solution and the enrichment of autologous blood. Manual use.

#### Indications

Extracorporeal withdrawal or transfer of blood and blood components and electrolyte solution for the preparation of blood and blood components used in the treatment of related diseases and conditions.

#### Contraindications

The product is not designed and approved for:

- Home use
- Use by unqualified, untrained personnel

The current Summary of Product Characteristics of the pharmaceutical manufacturers for the extracorporeal withdrawal, preparation (mixing, reconstitution) or transfer of blood/blood components/electrolyte solution with regard to incompatibilities (interaction of material/pharmaceuticals and interaction of different pharmaceuticals) and contraindications shall be observed.

#### Clinical benefits

The product in this scope enables the safe withdrawal, preparation (mixing, reconstitution) or transfer of blood and blood components. The clinical benefit of the product solely comes into effect when the prepared blood or blood component is administered. For administration, an appropriate medical device needs to be chosen.

#### General safety instructions

• The product is a medical device. • The product shall only be used by qualified medical and pharmaceutical personnel in hospitals, specialised practices and pharmacies with conventional environmental and hygiene conditions. • The product shall not be used on patients. • The instructions for use shall be observed before using the product. If applicable, further accompanying product information shall be observed before use. • If the sterile packaging is damaged or unintentionally opened before use, the product must be completely disposed of and replaced by a new one. • The product is intended for single use. Resterilisation as well as reuse are not permitted. Reuse can cause infections. • In the case of malfunctions (e.g. leakage, detached protective caps) or performance changes, the product shall not be used, and must be replaced. The defective product is to be provided to the manufacturer in compliance with measures for safe transport. • Due to the sterilisation method, the product can contain traces of ethylene oxide and its residual products. In case of known allergies to these substances, refrain from using the product. • The international and

national guidelines for the extracorporeal withdrawal, preparation (mixing, reconstitution) or transfer of blood and blood components shall be taken into account.

• Always observe the instructions for use from the pharmaceutical manufacturers regarding incompatibilities (interaction material/drugs as well as interaction of different drugs) and the current Summary of Product Characteristics. Material lists are available on request for all medical devices manufactured by CODAN. • The function of the product shall be checked by qualified medical/pharmaceutical personnel at regular intervals during use. • Generally accepted hygiene measures and working methods must be observed at all times. • The use of protective medical gloves is recommended to prevent any infections of users.

#### Notes on safe use of product

• The product is designed for manual use. • Prior to use, the product shall be checked for damage, tightness and any misconnections. • Prior to use, existing closure plugs and disconnectable connections must be checked for tightness and retightened if necessary. • Prior to and during use of the product, all connection ports shall be disinfected before and after each connection/disconnection. Isopropyl alcohol, iodine-containing alcohol or ethyl alcohol 70% or isopropyl alcohol 70% in combination with chlorhexidine gluconate 2% can be used for disinfection. Please observe that the exposure time has to amount to a minimum of 30 seconds. It is to be ensured that a connection is only made when completely dried. • If applicable, please ensure that the wall of the blood bag remains intact during piercing. • Excessive forces on Luer-Lock connections can lead to stress cracks.

Product with SWAN-LOCK®:

• The SWAN-LOCK® can be used for up to 140 activations during a permitted period of use. • Do not use cannulas or similar in conjunction with the SWAN-LOCK®. • The use of Luer-Lock connectors is recommended for a safe connection. • When connecting a male Luer slip cone, a simultaneous push and rotary movement is required in order to ensure a secure connection and a correct valve function. Do not use Luer connectors with excess length. These can impair or damage the function of the valve. • If possible, rinse the SWAN-LOCK® with a neutral solution after each application or withdrawal. Residues are rinsed out of the SWAN-LOCK®, preventing possible incompatibilities and securing the function of the valve against sticking due to residues. • The SWAN-LOCK® closes automatically when the syringe/line is disconnected. No sterile sealing cap is required.

#### Product compatibility

• The closure-piercing device complies with EN ISO 1135-4 and can be connected with all compatible infusion/blood containers. • All Luer and Luer-Lock connectors are compatible according to EN ISO 80369-1/ EN ISO 80369-7.

## Residual risks

No significant residual risks are known at this time.

## Period of use

The international and national guidelines for the extracorporeal withdrawal, preparation (mixing, reconstitution) or transfer of blood and blood components shall be taken into account.

The product is designed for a period of use of maximum 12 hours.

## Product disposal

The product shall be disposed of in accordance with the medical regulations for the disposal of contaminated medical waste. Protective medical gloves must be worn during disposal.

## Notification of serious incidents

All serious incidents occurring in connection with the product shall be reported to the manufacturer and the corresponding national reporting authority without delay. Defective products are to be provided to the manufacturer in compliance with measures for safe transport.

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