INSTRUCTIONS FOR USE

Cautery Tip Cleaner

Non-Sterile

Aspen Surgical Products Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316
Phone +1 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

EC REP

Emgro Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

For Country of Origin See Product Label

English

Intended Use

Intended to assist in cleaning and removing material from cautery tips during surgical procedures.

Intended User/Patient Target Groups

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Instructions for Use:

- Remove backing and apply to appropriate surface.
- Wipe cautery tip on pad to remove debris.
- Device is intended for limited contact duration, less than or equal to 16 hours.

Sterilization Instructions:

- Place device into a kit or sterilization tray.
- Using Ethylene Oxide sterilization method, validate according to ISO 11135, sterilize according to the following parameters:

Cycle Parameter	Specification	Duration
Preconditioning		
Temperature	110°F +/- 10°F	
Humidity	65% +/- 15% RH	10.35 Hrs - 72 Hrs
Chamber		
Vacuum	21" +/- 1.0 inHgA	
Humidification/Dwell	22.5" +/- 1.0 HgA	60 +30/-0 minutes
Gas Inject/Dwell	10.3" +/- 0.5 Hg rise	360 +5/-0 minutes
Jacket Temperature	135° +/- 5° F	
Aeration*		
Temperature	120° +/- 10° F	24 Hrs - 72 Hrs Max

*Proper aeration is essential to eliminating the risk of residuals to meet ISO 10993-7 requirements.

Warning

Reuse of device could result in infection/contamination and/or device failure which could lead to patient harm.

Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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AB'deki Kullanıcılar ve/veya Hastalara Yönelik Bildirim
Emergo Consulting (UK) Limited