


**LEGEND:**

	BREAK LINE
	CENTER LINE
	CUT LINE
	DIMENSION
	FOLD LINE
	1/8" INK BLEED
	PERFORATED LINE
	PLANT NOTES
	PRODUCT
	TEAR LINE
	TEXT
	TITLE BLOCK
	ZIPPER LINE

**NOTES:**

ARTWORK PRINT PROCESS: OFFSET
VARIABLE PRINT PROCESS: INKJET
VARIABLE PRINT INFORMATION CONTAINS REF. NUMBER, LOT NUMBER, EXPIRATION DATE, NEEDLE GAGE, NEEDLE LENGTH AND 2D BARCODE.
MUST RESPECT SIZE, LOCATION AND ORIENTATION OF VARIABLE AREAS.
CARTON IS VARNISHED EXCEPT WHERE NOTED WITH "NO VARNISH" HATCH.
*.125 MARGIN FROM ALL PRINT AREAS TO ALL EDGES, FOLDS AND SCORES UNLESS NOTED OTHERWISE.



<b>COVIDIEN</b> 15 Hampshire Street, Mansfield, MA 02048			
Job #:	MS2465	Business Unit:	Medical Supplies
Revision Date:	02/15/2018	Artwork/Part #:	DL29515138703
File Name:	DL29515138703_GNRC_INTR		
Template ID:	DL29515138702_GNRC_INTR_COVTMP		
Dimensions:	8.000" W x 2.000" H x 3.250" D	Outline Name:	5151387-U_REV0_PR
Print Process:	Offset	Package/Item Type:	INTR (Secondary) Carton
Designer(s):	K. Monteiro	Application/Version:	InDesign/CC
Associated to:	8881401031, 8881401056		
Comments/Printer Notes:	The use of the 2D Data Matrix will be permitted on the Cartons (display boxes) in Augusta, Norfolk and Deland. TLU unchanged from existing artwork due to space constraints and proximity to perforation. [2.15.18 KM - Updated symbols]		