

Instructions for Use for BARDIA™ Bed Bags

Contains or Presence of Phthalates: Di(2-ethylhexyl) phthalate (DEHP) is a plasticizer used in some polyvinyl chloride medical devices. DEHP has been shown to produce a range of adverse effects in experimental animals, notably liver toxicity and testicular atrophy. Although the toxic and carcinogenic effects of DEHP have been well established in experimental animals, the ability of this compound to produce adverse effects in humans is controversial. There is no evidence that neonates, infants, pregnant and breast feeding women exposed to DEHP experience any related adverse effects. However, a lack of evidence of causation between DEHP-PVC and any disease or adverse effect does not mean that there are no risks.

A space has been provided for you to write the date on which the urine bag should be changed. Please follow local policy regarding how long the bag can remain in place.

This is a single use device. Do not resterilize any portion of this device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.

BARD



Manufacturer:
C. R. Bard, Inc.
Covington, GA 30014 USA
1-800-526-4455
www.bardmedical.com

EC REP

Bard Limited
Forest House
Brighton Road
Crawley, West Sussex UK
RH11 9BP
+44 1293 527 888

CE
0086

PK7646642 10/2017